

## Bacterial Colony Counter


### Application:

Widely used for bacilli inspection of food, drink, medicine, biological, make-up, health products, Water for drink, domestic sewage, and industrial water.

The required equipment of all kinds' health station, environmental monitoring centre, food hygiene, testing centre, hospital, Institute of biological products, laboratory for the control of drugs, Commodity Inspection bureau, the made-up factory, the laboratory of university and scientific research institutions.

### Features:

- \* LED Digital display
- \* Made of counter, pen, counting cell. The counter is made of CMOS integrated electric, LED Nixie tube.
- \* The word high is 13mm , clear and bright; The pen is sensitive and accurate
- \* There are two kinds of light source, adjustable.



BC-6

BC-6A

### Technical Parameters:

Model	Specification
BC-6	<p>Automatic, Optical amplification 8 x, reflection and transmission illumination</p> <p>Counter range: 1~999</p> <p>Suitable petri dishes: ø50mm, ø90mm</p> <p>External size: 310*210*80mm</p> <p>Power supply: AC110/220V±10%, 50/60Hz</p> <p>Package size: 510*410*280mm</p> <p>Gross weight: 2.5kg</p>
BC-6A	<p>Automatic, Optical amplification 5 x 10 x, reflection and transmission illumination</p> <p>Counter range: 1~999</p> <p>Suitable petri dishes: ø50mm, ø90mm</p> <p>External size: 310*210*260mm</p> <p>Power supply: AC110/220V±10%, 50/60Hz</p> <p>Package size: 510*410*460mm</p> <p>Gross weight: 3.8kg</p>